

## BALL-HEAD

The Champions® Implant Ball-Head is a screw-shaped endosseous one-piece implant. The implant has diameters of 2.5/3.0/4.0 mm. It is available in lengths of 8/10/12/14/16/18 mm. The implant covers a wide range of indications.

The implant is suitable for the following indications:

Insertion in the long-term healed jaw bone (delayed implantation)


Delayed insertion (6-8 weeks after the tooth extraction)

Immediate implantation (implantation directly after the tooth extraction, under certain conditions)

The Champions® Implant Ball-Head can be placed in bone in the maxilla and mandible (D1-D4) to support a removable prosthesis. There should be at least 3 implants/teeth per quadrant in the maxilla and at least 2 implants/teeth per quadrant in the mandible.

Please note the restrictions of indications in the Instructions for Use.

Material: Titanium grade 4



Diameter	Length	Height	Art. No.	Type
2.5 mm	10 mm	6.3 mm	1900	●
	12 mm	6.3 mm	2015	
	14 mm	6.3 mm	2035	
	16 mm	6.3 mm	2055	
3.0 mm	8 mm	6.8 mm	1890	
	10 mm	6.8 mm	1905	
	12 mm	6.8 mm	2016	
	14 mm	6.8 mm	2040	
	16 mm	6.8 mm	2060	
	18 mm	6.8 mm	2080	
4.0 mm	8 mm	7.1 mm	1895	
	10 mm	7.1 mm	2005	
	12 mm	7.1 mm	2025	
	14 mm	7.1 mm	2050	
	16 mm	7.1 mm	2070	