

IMPLANTS

SQUARE "CLASSIC"

The Champions® Implant Square Classic is a screw-shaped endosseous one-piece implant. The implant has a diameter of 3.5 mm. It is available in lengths of 8/10/12/14/16 mm. The implant covers a wide range of indications.

The implant is suitable for the following indications:

Insertion in the long-term healed jaw bone (delayed implantation)


Delayed insertion (6-8 weeks after the tooth extraction)

Immediate implantation (implantation directly after the tooth extraction, under certain conditions)

The Champions® Implant Square Classic can be placed in bone in the mandible (D1).

Please note the restrictions of indications in the Instructions for Use.

Material: Titanium grade 4



Diameter	Length	Height	Art. No.	Type
3.5 mm	8 mm	8.1 mm	1005	●
	10 mm	8.1 mm	1025	
	12 mm	8.1 mm	1055	
	14 mm	8.1 mm	1095	
	16 mm	8.1 mm	1135	